


**UNIQUE  
CHALLENGES.  
INDIVIDUALIZED  
SOLUTIONS.**

PORTFOLIO OF SOLUTIONS  
FOR COMPLEX PCI



# DRIVING SUCCESSFUL OUTCOMES FOR HIGH-RISK PCI PATIENTS

Every complex case is uniquely challenging. That's why we focus on providing the therapies, tools, and specialized knowledge that enable you to take an individualized approach and optimize long-term durable outcomes.

## ENSURING MAXIMUM PREPARATION

We bring notable clinical acumen and product expertise to support optimal clinical preparation in advance of treatment.

## INDIVIDUALIZED INNOVATION & GUIDANCE

We invest in novel tools and deliver lesion-level specialized guidance to meet individualized needs in your most challenging cases.

## IMPROVING LONG-TERM OUTCOMES


We advance the latest evidence-based solutions to improve clinical outcomes and economic value for the most challenging interventional procedures.



# TOOLS TO TAKE ON THE COMPLEX

We equip you with a growing portfolio of therapies and tools for complex cases — from an innovative non-compliant scoring balloon, catheters for tracking tortuous anatomy, and ultra-low profile microcatheters for challenging lesions to the proven dual-action efficacy of our orbital atherectomy system.

## SCOREFLEX® NC SCORING PTCA CATHETER\*




scoreflex® NC  
Scoring PTCA Catheter

## FOCUSED-FORCE ANGIOPLASTY

Scoreflex® NC is a focused-force dilatation balloon for non-compliant plaque modification and lesion preparation. It uses a dual-wire system, which creates a focal stress pattern to facilitate safe and controlled plaque modification at lower resolution pressure.

- **First non-compliant coronary scoring balloon**
- **Highest-rated burst pressure**
- **Low crossing profile**
- **Indicated for ISR**

## TELEPORT® MICROCATHETER\*



Teleport®  
Microcatheter


## ULTIMATE CONTROL AND TRACKABILITY

- **Ultra low profile and pushability** for navigation inside challenging lesions and micro channels
- **HYbrid BRAiding + Coil (HYBRACOIL) construction** for guidewire support and crossing
- **DUal-Layer RAdiopaque (DURA) tapered tip** for robustness and visibility
- **Indicated for coronary and peripheral intervention**

\*Manufactured by:

 **OrbusNeich™**  
Pioneers in life-changing technologies

## SAPPHIRE® II PRO DILATATION CATHETER\*




**SAPPHIRE® II PRO**  
Balloon Dilatation Catheter

## CROSS WITH CONFIDENCE

The Sapphire® II PRO is specifically engineered for crossing difficult lesions and tracking tortuous anatomy.


- **Featuring the first and only 1.0 mm balloon with FDA clearance\*\***
- **Excellent pushability**
- **Robust material**
- **Ultimate crossability**

**Sub-zero tip technology for smooth lesion entry**



The well-balanced, sub-zero tapered tip has an ultra-low profile and provides smooth entry through the tightest lesions.

## SAPPHIRE® NC 24 NON-COMPLIANT PTCA CATHETER\*



**SAPPHIRE® NC 24**  
Coronary Dilatation Catheter

## CONTROLLED COMPLIANCE

- **Highest RBP on the US market at 24 atm\*\***
- **High-pressure tolerance and recrossability**
- **Controlled balloon compliance** for the most accurate sizing and stent optimization

\*Manufactured by:

 **OrbusNeich™**  
Pioneers in life-changing technologies

\*\* at time of publication

# SEVERE CALCIUM ONE SOLUTION

The Diamondback 360® Orbital Atherectomy System is your single solution for severely calcified coronary artery disease, with proven success in challenging anatomies and lesions.

**VIPERWIRE ADVANCE® NITINOL**  
coronary guidewire with flex tip<sup>1</sup>



**ONE-TOUCH START BUTTON**  
makes device power up effortless

**GLIDEASSIST®**

Easier tracking and removal and smoother repositioning of the device – especially in challenging anatomies†



**SINGLE 1.25 MM CROWN**  
treats vessels 2.5 to 4.0

**ELECTRIC-POWERED HANDLE**  
allows two-minute set up and provides efficient torque transfer to the shaft and crown

**TWO CONVENIENT SPEED CONTROLS**  
allow for quick speed adjustments within the sterile field

## DUAL ACTION

Uniquely designed to enable simultaneous modification of both intimal and medial calcium for optimal stent delivery, expansion, and apposition in severely calcified lesions.

## VERSATILE

Treat a broad range of challenging cases, including the most severely calcified lesions, with under 2-minute setup and predictable procedure times.

## PROVEN

Extensively studied, with over 145,000 patients treated, demonstrated to perform effectively in the treatment of severely calcified lesions.

† as compared to use of Diamondback without GlideAssist



TO SEE ACCESSORIES,  
SCAN THE QR CODE



# PRODUCT ORDERING INFORMATION

**DIAMONDBACK 360°**  
CORONARY ORBITAL ATHERECTOMY SYSTEM

## DIAMONDBACK 360° ORBITAL ATHERECTOMY SYSTEM

DIAMONDBACK 360° ORBITAL ATHERECTOMY DEVICE			
MODEL NUMBER	CROWN SIZE	SHAFT LENGTH	QUANTITY
DBEC-125	1.25 mm Classic	135 cm	1 each

VIPERWIRE ADVANCE® CORONARY GUIDE WIRES			
MODEL NUMBER	SIZE	WIRE LENGTH	QUANTITY
GWC-12325LG-FLP	.012"/0.014" Tip	325 cm	5 per box
GWC-12325LG-FT	.012"/0.014" Flex Tip	325 cm	5 per box


VIPERSLIDE® LUBRICANT		
MODEL NUMBER	DESCRIPTION	QUANTITY
VPR-SLD2	100 mL Package	10 bags per box

OAS PUMP		
MODEL NUMBER	DESCRIPTION	QUANTITY
SIP-3000	OAS Pump	1 each

**SAPPHIRE® NC24**  
Coronary Dilatation Catheter

## SAPPHIRE® NC24 NON-COMPLIANT PTCA CATHETER\*

BALLOON DIAMETER (mm)	BALLOON WORKING LENGTH (mm)						
	8	10	12	15	18	22	26
2.00	492000812	492001012	492001212	492001512	492001812	492002212	492002612
2.25	492250812	492251012	492251212	492251512	492251812	492252212	492252612
2.50	492500812	492501012	492501212	492501512	492501812	492502212	492502612
2.75	492750812	492751012	492751212	492751512	492751812	492752212	492752612
3.00	493000812	493001012	493001212	493001512	493001812	493002212	493002612
3.25	493250812	493251012	493251212	493251512	493251812	493252212	493252612
3.50	493500812	493501012	493501212	493501512	493501812	493502212	493502612
3.75	493750812	493751012	493751212	493751512	493751812	493752212	493752612
4.00	494000812	494001012	494001212	494001512	494001812	494002212	494002612
4.50	494500812	494501012	494501212	494501512	494501812	494502212	494502612
5.00	495000812	495001012	495001212	495001512	495001812	495002212	495002612



BALLOON DIAMETER (mm)		BALLOON WORKING LENGTH (mm)						
		5	8	10	12	15	20	30
OTW	1.00	501000512	501000812	501001012	–	501001512	–	–
Rx	1.00	210-053-5UU	210-083-5UU	210-103-5UU	–	210-153-5UU	–	–
	1.25	212-053-5UU	212-083-5UU	212-103-5UU	–	212-153-5UU	–	–
	1.50	–	–	215-103-5UU	215-123-5UU	215-153-5UU	215-203-5UU	–
	1.75	–	–	217-103-5UU	–	217-153-5UU	217-203-5UU	–
	2.00	–	–	220-103-5UU	220-123-5UU	220-153-5UU	220-203-5UU	–
	2.25	–	–	222-103-5UU	–	222-153-5UU	222-203-5UU	–
	2.50	–	–	225-103-5UU	225-123-5UU	225-153-5UU	225-203-5UU	225-303-5UU
	2.75	–	–	227-103-5UU	–	227-153-5UU	227-203-5UU	–
	3.00	–	–	230-103-5UU	230-123-5UU	230-153-5UU	230-203-5UU	230-303-5UU
	3.25	–	–	232-103-5UU	–	232-153-5UU	232-203-5UU	–
	3.50	–	–	235-103-5UU	–	235-153-5UU	235-203-5UU	235-303-5UU
	4.00	–	–	240-103-5UU	–	240-153-5UU	240-203-5UU	–

BALLOON DIAMETER (mm)	BALLOON WORKING LENGTH (mm)		
	10	15	20
1.75	617-104-1U	617-154-1U	617-204-1U
2.00	620-104-1U	620-154-1U	620-204-1U
2.25	622-104-1U	622-154-1U	622-204-1U
2.50	625-104-1U	625-154-1U	625-204-1U
2.75	627-104-1U	627-154-1U	627-204-1U
3.00	630-104-1U	630-154-1U	630-204-1U
3.50	635-104-1U	635-154-1U	635-204-1U
4.00	640-104-1U	640-154-1U	640-204-1U

MODEL NUMBER	PROFILE	LENGTH
220-13-1000U	2.0 Fr	135 cm
220-15-1000U	2.0 Fr	150 cm
221-13-1000U	2.1 Fr	135 cm
221-15-1000U	2.1 Fr	150 cm

1. CSI data on file: based on cadaver atherosclerotic lesions, porcine coronary lesions and graphite block test models.

**Indication:** The DIAMONDBACK 360® Coronary Orbital Atherectomy Systems (OAS) are percutaneous orbital atherectomy system indicated to facilitate stent delivery in patients with coronary artery disease (CAD) who are acceptable candidates for PTCA or stenting due to de novo, severely calcified coronary artery lesions. **Contraindications:** The OAS are contraindicated when the VIPERWIRE® guide wire cannot pass across the coronary lesion or the target lesion is within a bypass graft or stent. The OAS are contraindicated when the patient is not an appropriate candidate for bypass surgery, angioplasty, or atherectomy therapy, or has angiographic evidence of thrombus, or has only one open vessel, or has angiographic evidence of significant dissection at the treatment site and for women who are pregnant or children. **Warnings/Precautions:** Performing treatment in excessively tortuous vessels or bifurcations may result in vessel damage; The OAS were only evaluated in severely calcified lesions, A temporary pacing lead may be necessary when treating lesions in the right coronary and circumflex arteries; On-site surgical back-up should be included as a clinical consideration; Use in patients with an ejection fraction (EF) of less than 25% has not been evaluated. See the Instructions for Use before performing DIAMONDBACK 360 coronary orbital atherectomy procedures for detailed information regarding the procedure, indications, contraindications, warnings, precautions, and potential adverse events. For further information call CSI at 1-877-274-0901 and/or consult CSI's website at [www.csi360.com](http://www.csi360.com). **Caution:** Federal law (USA) restricts this device to sale by or on the order of a physician.

**Indication:** The Scoreflex NC Scoring PTCA Catheter is indicated for balloon dilatation of a de novo stenotic portion of a coronary artery in patients evidencing coronary ischemia for the purpose of improving myocardial perfusion, in-stent restenosis. **Contraindications:** The use of the Scoreflex NC Scoring PTCA Catheter is contraindicated in the following patient types: Patients with an unprotected left main coronary artery; Patients with coronary artery spasm in the absence of a significant stenosis. **Warnings:** When using this type of device, the following warnings should be observed: To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis. When the catheter is exposed to the vascular system, it should be manipulated while under high-quality fluoroscopic observation. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding. Applying excessive force to the catheter can result in tip or catheter breakage, catheter kink, or balloon separation. Do not twist the catheter shaft in excess of 180 degrees when the tip is constrained. Balloon pressure should not exceed the rated burst pressure (RBP) indicated on the package. **Precautions:** Never advance the Scoreflex NC Scoring PTCA Catheter without the guidewire extending from the tip. Do not use oil-based contrast medium, organic solvents or alcohols; there is a possibility of catheter leak, damage, or lubrication loss. The balloon deflation time has been established as 15 seconds based on in vitro bench testing results. Do not reinsert the PTCA catheter into the coil dispenser after procedural use. Discard all disposable devices used during this procedure per local requirements for medical device waste disposal. **Caution:** Federal law (USA) restricts this device to the sale by or on the order of a physician.

CSI, Diamondback 360, GlideAssist and ViperWire Advance are registered trademarks of Cardiovascular Systems, Inc. © 2023 Cardiovascular Systems, Inc. EN-7378.A 1222 Sapphire, Scoreflex and Teleport are registered trademarks of OrbusNeich Medical Group Holdings Ltd or its affiliates. OrbusNeich is a trademark of OrbusNeich Medical Group Holdings Ltd or its affiliates. © 2023 OrbusNeich Medical Group Holdings Ltd or its affiliates. All Rights Reserved. G-70-2125 Rev 01. All other trademarks cited herein are the trademarks of their respective owners.

**Distributed by:**



1225 Old Hwy 8 NW  
St. Paul, MN 55112

T: 651.259.1600  
877.274.0901

F: 612.677.3355  
W: [csi360.com](http://csi360.com)

**\*Manufactured by:**



1 Jinkui Road,  
Fuitian Free Trade Zone,  
Shenzhen, 518038, China

T: +86.755.8358.0181  
F: +86.755.8358.0169  
[www.OrbusNeich.com](http://www.OrbusNeich.com)